




Tiras de temperatura irreversibles de 8 pasos testoterm[®] **Testo**

Tiras de temperatura irreversibles de 8 pasos testoterm[®]

Tiras para 8 temperaturas. 5 tipos cubren el rango de medición de 37 a 260 °C, cada tira muestra 8 temperaturas diferentes. Cada tira mide 50 x 18 mm.



Code	Description	Packaging
LLG09246002	Temperature recording strips, 51 x 18/10	10 pz.
LLG09246003	Temperature recording strips, 51 x 1/10	10 pz.
LLG09246004	Temperature recording strips, 51 x 1/10	10 pz.
LLG09246005	Temperature recording strips, 51 x 1/10	10 pz.
LLG09246006	Temperature recording strips, 51 x 1/10	10 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

