



Tiras de temperatura irreversibles de 5 pasos Thermax[®]

Tiras de temperatura irreversibles de 5 pasos Thermax[®]

Tiras de medición de temperatura irreversibles de 5 pasos que cubren el rango de temperatura de 29 a 290°C. Cada tira mide 35 x 15 mm.



Code	Description	Packaging
LLG09246310	Temperature recording strips, 35 x 1/10	10 pz
LLG09246311	Temperature recording strips, 35 x 1/10	10 pz
LLG09246312	Temperature recording strips, 35 x 1/10	10 pz
LLG09246313	Temperature recording strips, 35 x 1/10	10 pz
LLG09246314	Temperature recording strips, 35 x 1/10	10 pz
LLG09246315	Temperature recording strips, 35 x 1/10	10 pz
LLG09246316	Temperature recording strips, 35 x 1/10	10 pz
LLG09246317	Temperature recording strips, 35 x 1/10	10 pz
LLG09246318	Temperature recording strips, 35 x 1/10	10 pz
LLG09246319	Temperature recording strips, 35 x 1/10	10 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

