

**Embudos, DURAN<sup>®</sup>****DURAN Group****Embudos, DURAN<sup>®</sup>**

DIN ISO 4798. Con vástago corto. Tipo de vidrio I/vidrio neutro según USP, EP y JP. Autoclavables.



Code	Description	Packaging
LLG09251023	Funnels 120mm DURANDURAN glass	1 pz
LLG09251036	Funnels 180mm DURAN	1 pz.
LLG09251123	Funnels,DURAN,short stem,diam. 35 mm	1 pz
LLG09251128	Funnels,DURAN,short stem,diam. 45 mm	1 pz
LLG09251133	Funnels,DURAN,short stem,diam. 55 mm	1 pz
LLG09251138	Funnels,DURAN,short stem,diam. 70 mm	1 pz
LLG09251141	Funnels,DURAN,short stem,diam. 80 mm	1 pz
LLG09251146	Funnels,DURAN,short stem,diam. 100 mm	1 pz
LLG09251157	Funnels,DURAN,short stem,diam. 150 mm	1 pz
LLG09251161	Funnels,DURAN,short stem,diam. 200 mm	1 pz
LLG09251169	Funnels 300mm DURAN	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

