



**Embudos, acero inoxidable Remanit® 4301**

RSG Rostfrei-  
Schneidwerkzeuge

**Embudos, acero inoxidable Remanit® 4301**  
Con canal de aireación y asa.



Code	Description	Packaging
LLG04661801	Funnel 100 mm stainless steel,	1 pz
LLG04661834	Funnel 250 mm stainless steel	1 pz
LLG04661835	Funnel stainlesssteel 200mm handle + air	1 pz.
LLG04661836	Funnel stainlesssteel 150mm handle + air	1 pz.
LLG04664420	Funnel 50 mm, stainless steel, w/o handl	1 pezzo



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

