



Embudos, acero inoxidable Remanit® 4301

RSG Rostfrei-
Schneidwerkzeuge

Embudos, acero inoxidable Remanit® 4301
Con canal de aireación y asa.



Code	Description	Packaging
LLG04661801	Funnel 100 mm stainless steel,	1 pz
LLG04661834	Funnel 250 mm stainless steel	1 pz
LLG04661835	Funnel stainlesssteel 200mm handle + air	1 pz.
LLG04661836	Funnel stainlesssteel 150mm handle + air	1 pz.
LLG04664420	Funnel 50 mm, stainless steel, w/o handl	1 pezzo



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

