

**Embudos para análisis, DURAN[®]**

Para filtraciones rápidas. Tipo de vidrio I/vidrio neutro según USP, EP y JP.

Autoclavables.



Code	Description	Packaging
LLG09251237	Analytical funnels,DURAN,diam. 65 mm	1 pz
LLG09251241	Analytical funnels,DURAN,diam. 80 mm	1 pz
LLG09251248	Analytical funnels,DURAN,diam. 110 mm	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

