

**Embudos, PP****VITLAB****Embudos, PP**

Transparentes. Flujo rápido debido a un ángulo inclinado de 60°. Mango práctico con asa para colgar. Aptos para alimentos conforme a la normativa CE nr. 10/2011



| Code        | Description                            | Packaging |
|-------------|----------------------------------------|-----------|
| LLG09251602 | Funnels,PP,diam. 30 mm                 | 1 pz.     |
| LLG09251605 | Funnels,PP,diam. 30 mm,stem diam. 5 mm | 1 pz.     |
| LLG09251610 | Funnels,PP,diam. 40 mm                 | 1 pz.     |
| LLG09251615 | Funnels,PP,diam. 50 mm                 | 1 pz.     |
| LLG09251620 | Funnels,PP,diam. 75 mm                 | 1 pz.     |
| LLG09251625 | Funnels,PP,diam. 100 mm                | 1 pz.     |
| LLG09251630 | Funnels,PP,diam. 120 mm                | 1 pz.     |
| LLG09251635 | Funnels,PP,diam. 150 mm                | 1 pz.     |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6ª Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

