

**Crown stopper, mechanical****Hanhart****Crown stopper, mechanical**

Impact resistant, dustproof and water resistant. With diamond-turned metal case.  
55mm diameter. Built-in strap ring.

**CARLO ERBA**  
REAGENTS

Code	Description	Packaging
LLG09262005	Crown stopper, mechanical, Type Crown stopper 1 / 5 s, Measuring capacity 00:30:00 hr: min: sec	1 pz
LLG09262010	Crown stopper, mechanical, Type Crown stopper 1 / 10 s, Measuring capacity 00:15:00 hr: min: sec	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

