



Vidrios de reloj, DURAN[®]

DURAN Group

Vidrios de reloj, DURAN[®]

- DIN 12341
- Reborde fundido
- Esterilizable en autoclave



Code	Description	Packaging
LLG09263424	Watch glasses,DURAN,diam. 40 mm	1 pz
LLG09263432	Watch glasses,DURAN,diam. 50 mm	1 pz
LLG09263434	Watch glasses,DURAN,diam. 60 mm	1 pz
LLG09263441	Watch glasses,DURAN,diam. 80 mm 1 pz	1 pz
LLG09263446	Watch glasses,DURAN,diam. 100 mm	1 pz
LLG09263452	Watch glasses,DURAN,diam. 125 mm	1 pz
LLG09263457	Watch glasses,DURAN,diam. 150 mm	1 pz
LLG09263461	Watch glasses,DURAN,diam. 200 mm	1 pz
LLG09263466	Watch glasses,DURAN,diam. 250 mm	1 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5ª Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

