



**Vidrios de reloj, DURAN<sup>®</sup>**

**DURAN Group**

**Vidrios de reloj, DURAN<sup>®</sup>**

-DIN 12341

-Reborde fundido

-Esterilizable en autoclave



Code	Description	Packaging
LLG09263424	Watch glasses,DURAN,diam. 40 mm	1 pz
LLG09263432	Watch glasses,DURAN,diam. 50 mm	1 pz
LLG09263434	Watch glasses,DURAN,diam. 60 mm	1 pz
LLG09263441	Watch glasses,DURAN,diam. 80 mm 1 pz	1 pz
LLG09263446	Watch glasses,DURAN,diam. 100 mm	1 pz
LLG09263452	Watch glasses,DURAN,diam. 125 mm	1 pz
LLG09263457	Watch glasses,DURAN,diam. 150 mm	1 pz
LLG09263461	Watch glasses,DURAN,diam. 200 mm	1 pz
LLG09263466	Watch glasses,DURAN,diam. 250 mm	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

