



Viscosímetro de inmersión con boquilla fija

Viscosímetro de inmersión con boquilla fija
Con boquilla fija de acero inoxidable o aluminio.



Code	Description	Packaging
LLG06080098	Dipping flow cups TA, aluminum with fix	1 pz.
LLG06225937	Viscosity cup Type TAR, stainless steel	1 pz.
LLG06264682	Viscosity cup Type TAR, stainless steel	1 pz.
LLG06302328	Viscosity cup Type TI, aluminium with f	1 pz.
LLG06314566	Viscosity cup Type TI, aluminium with f	1 pz.
LLG09149894	Dipping flow cups TA aluminum with 4 mm	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

