



## Viscosímetro de inmersión con boquilla fija

**Viscosímetro de inmersión con boquilla fija**  
Con boquilla fija de acero inoxidable o aluminio.



Code	Description	Packaging
LLG06080098	Dipping flow cups TA, aluminum with fix	1 pz.
LLG06225937	Viscosity cup Type TAR, stainless steel	1 pz.
LLG06264682	Viscosity cup Type TAR, stainless steel	1 pz.
LLG06302328	Viscosity cup Type TI, aluminium with f	1 pz.
LLG06314566	Viscosity cup Type TI, aluminium with f	1 pz.
LLG09149894	Dipping flow cups TA aluminum with 4 mm	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

