



Wide-mouth bottles, clear and amber glass, PTFE-lined screw caps

behr

Wide-mouth bottles, clear and amber glass, PTFE-lined screw caps

Clear and amber glass. For storage and transport. The PTFE-lined screw caps make them ideal for taking samples.



Code	Description	Packaging
LLG04658141	Wide neck bottle RB 30 GT, 30ml amber/1	1 pz
LLG04658142	Wide neck bottle RB 50 GT, 50ml amber/1	1 pz
LLG04658143	Wide neck bottle RB 100 GT, 100ml amb/1	1 pz
LLG04658144	Wide neck bottle RK 30 GT, 30ml clear/1	1 pz
LLG04658145	Wide neck bottle RK 50 GT, 50ml clear/1	1 pz
LLG09072325	Wide mouth bottle 250ml, amberwith PTFE	1 pz
LLG09072326	Wide mouth bottles 500ml amber + closure	1 pz
LLG09072327	Wide-mouth bottles, clear and amber glass, PTFE-lined screw caps, Type RB 1000 GT, Colour amber, 1000 ml, Thread diam. 68 mm	1 pz
LLG09103200	Wide-neck bottle RK 250 GT250 ml, threa	1 pz
LLG09103201	Wide-neck bottle 500 ml clear + closure	1 pz
LLG09103202	Wide-mouth bottles, clear and amber glass, PTFE-lined screw caps, Type RK 1000 GT, Colour clear, 1000 ml, Thread diam. 68 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

