



Recipiente para almacenar pipetas

Kartell

Recipiente para almacenar pipetas

Con tapa. Ideal para proteger del polvo pipetas de hasta 400 mm de longitud.



Code	Description	Packaging
LLG09273900	Pipet storage container,PP,130 x 430 mm	1 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

