



**Recipiente para almacenar pipetas**

**Kartell**

**Recipiente para almacenar pipetas**

Con tapa. Ideal para proteger del polvo pipetas de hasta 400 mm de longitud.



Code	Description	Packaging
LLG09273900	Pipet storage container,PP,130 x 430 mm	1 pz



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

