

**Tarros de boca ancha, PET, con tapa de rosca, PP****Behr****Tarros de boca ancha, PET, con tapa de rosca, PP**

Transparentes con tapa de polipropileno blanco. Aptos para uso alimentario.  
Redondos, excepto 1000 ml y 2000 ml que son cuadrados.



Code	Description	Packaging
LLG06242746	Jar, PET, wide-mouth PWG 1001 clear tra	10 pz.
LLG09103136	PET wide mouth bottlespack of 10, 10/10	10 pz.
LLG09103137	PET wide mouth bottlespack of 10, 25/10	10 pz.
LLG09103138	PET wide mouth bottles 500mlpack of /10	10 pz.
LLG09103140	PET wide mouth bottles 2000mlpack of/10	10 pz.
LLG09103148	PET wide mouth bottles 2500mlpack of/10	10 pz.
LLG09103149	PET wide mouth bottles 3000mlpack of/10	10 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

