



Jarras graduadas, con asa, de acero inoxidable, alta resistencia

Bochem

Jarras graduadas, con asa, de acero inoxidable, alta resistencia

De acero 18/10. Graduadas 100 ml. Con pico y asa.



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09275850 | Measuring beaker cap. 1 ltr.18/10 steel | 1 pz. |
| LLG09275852 | Measuring beaker cap. 2 ltr.18/10 steel | 1 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

