



Accesorios para Titronic[®] 300 y TitroLine[®] 5000



Code	Description	Packaging
LLG09863188	Magnetic stirrer TM 50 for glass vessel	1 pz.
LLG09863189	USB-HUB TZ 3830, to connect more device	1 pz.
LLG09863190	Reagent bottle 1000ml, TZ 3803 amber gl	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

