


**Probetas de PP, clase B, graduación en relieve****Kartell****Probetas de PP, clase B, graduación en relieve**

Forma alta, con pico, alta claridad. Conforme a ISO 6706 - 1981 (E). Calibradas para contener (TC, In). Aceptan la exposición a temperaturas de hasta 80 °C, que no les provoca el sobrepaso permanente de los límites de tolerancia. Se recomienda una temperatura de limpieza inferior a 60 °C para conservar las marcas e inscripciones. Aptas para alimentos, esterilizables en autoclave.



Code	Description	Packaging
LLG09274750	Measuring cylinder 10ml, tall formblue	1 pz.
LLG09274751	Measuring cylinder 25ml, tall formblue	1 pz.
LLG09274752	Measuring cylinder 50ml, tall formblue	1 pz.
LLG09274753	Measuring cylinder 100ml, tall formblue	1 pz.
LLG09274755	Measuring cylinder 500ml, tall formblue	1 pz.
LLG09274756	Measuring cylinder 1000ml, tall formblu	1 pz.
LLG09274757	Measuring cylinder 2000ml, tall formcbl	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

