




Tapones de aluminio para encapsulado ND8 LLG, ensamblados

LLG Labware

Tapones de aluminio para encapsulado ND8 LLG, ensamblados



Code	Description	Packaging
LLG04008198	LLG-Aluminium crimp cap N 8, silver/100	100 pz
LLG04008200	LLG-Aluminium crimp cap N 8, silver/100	100 pz
LLG06266865	LLG-Alu-Vials closure N8, 45° shore/100	100 pz.
LLG09003443	LLG-Vials closure N8 TB/oA, 1,3mmpa/100	100 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

