




Tapones de aluminio para encapsulado ND8 LLG, ensamblados

LLG Labware

Tapones de aluminio para encapsulado ND8 LLG, ensamblados



Code	Description	Packaging
LLG04008198	LLG-Aluminium crimp cap N 8, silver/100	100 pz
LLG04008200	LLG-Aluminium crimp cap N 8, silver/100	100 pz
LLG06266865	LLG-Alu-Vials closure N8, 45° shore/100	100 pz.
LLG09003443	LLG-Vials closure N8 TB/oA, 1,3mmpa/100	100 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

