



Cierres de rosca de PP ND8, ensamblados, y tapones de rosca vacíos ND8 LLG

LLG Labware

Cierres de rosca de PP ND8, ensamblados, y tapones de rosca vacíos ND8  
LLG



Code	Description	Packaging
LLG04008207	LLG-Screw cap N 8, PP, black, cente/100	100 pz
LLG04008208	LLG-Screw cap N 8, PP, black, close/100	100 pz
LLG06232178	LLG-screw-caps N 8, blue,with head /100	100 pz
LLG06266866	LLG-Screw cap 8 mm, PP, blackwith h/100	100 pz.
LLG06266867	LLG-Screw cap 8 mm, PP, blackclosed/100	100 pz.
LLG07060421	LLG-Screw cap N 8, PP, black, cente/100	100 pz
LLG07075960	LLG-Screw cap N 8, PP, black, closed to	100 pz
LLG09003484	LLG-Sscrew cap N 8, PP, black, cent/100	100 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

