



Cierres de aluminio para encapsulado ND11 LLG, ensamblados

LLG Labware

Cierres de aluminio para encapsulado ND11 LLG, ensamblados



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG04001522 | LLG-Aluminium crimp caps N 11 TB/oAgree | 100 pz    |
| LLG04001559 | LLG-Aluminium crimp cap N 11, silve/100 | 100 pz    |
| LLG06291637 | LLG-Aluminium crimp cap N 11, silverred | 100 pz    |
| LLG06900233 | LLG-Aluminium caps N 11 TB/oAwith s/100 | 100 pz    |
| LLG07050759 | LLG-Aluminium crimp cap N 11, silve/100 | 100 pz    |
| LLG07060469 | LLG-Aluminium crimp cap N 11, silver, c | 100 pz    |
| LLG07300348 | LLG-Alu flange caps N 11 TB/oAwith /100 | 100 pz    |
| LLG09003441 | LLG-Aluminium crimp cap N 11, silver, c | 100 pz    |
| LLG09003446 | LLG-Aluminium crimp cap N 11, silve/100 | 100 pz    |
| LLG09003459 | Alu rim caps N 11, goldhole cap, Bu/100 | 100 pz.   |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5ª Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

