



Cierres de aluminio para encapsulado ND11 LLG, ensamblados

LLG Labware

Cierres de aluminio para encapsulado ND11 LLG, ensamblados



Code	Description	Packaging
LLG04001522	LLG-Aluminium crimp caps N 11 TB/oAgree	100 pz
LLG04001559	LLG-Aluminium crimp cap N 11, silve/100	100 pz
LLG06291637	LLG-Aluminium crimp cap N 11, silverred	100 pz
LLG06900233	LLG-Aluminium caps N 11 TB/oAwith s/100	100 pz
LLG07050759	LLG-Aluminium crimp cap N 11, silve/100	100 pz
LLG07060469	LLG-Aluminium crimp cap N 11, silver, c	100 pz
LLG07300348	LLG-Alu flange caps N 11 TB/oAwith /100	100 pz
LLG09003441	LLG-Aluminium crimp cap N 11, silver, c	100 pz
LLG09003446	LLG-Aluminium crimp cap N 11, silve/100	100 pz
LLG09003459	Alu rim caps N 11, goldhole cap, Bu/100	100 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

