



Cierres de PE con aro de presión ND11 LLG, ensamblados

LLG Labware

Cierres de PE con aro de presión ND11 LLG, ensamblados



Code	Description	Packaging
LLG04001544	LLG-Snap ring cap N 11,PE, transpar/100	100 pz
LLG04008256	LLG-Snap ring cap N 11, PE, transparent	100 pz
LLG04008257	LLG-Snap ring cap N 11, PE, blue, c/100	100 pz
LLG04008258	LLG-Snap ring cap N 11PE, blue, center	100 pz
LLG04008259	LLG-Snap ring cap N 11, PE, bluecenter	100 pz
LLG04008261	LLG-Snap ring cap N 11, PEtransparent,c	100 pz
LLG06073555	PE Snap Ring Seals ND11,ready assembled	100 pz
LLG06242212	LLG-Snap ring cap N 11blue, hole, s/100	100 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

