



Cierres de PE con aro de presión ND11 LLG, ensamblados

LLG Labware

Cierres de PE con aro de presión ND11 LLG, ensamblados



Code	Description	Packaging
LLG04001544	LLG-Snap ring cap N 11,PE, transpar/100	100 pz
LLG04008256	LLG-Snap ring cap N 11, PE, transparent	100 pz
LLG04008257	LLG-Snap ring cap N 11, PE, blue, c/100	100 pz
LLG04008258	LLG-Snap ring cap N 11PE, blue, center	100 pz
LLG04008259	LLG-Snap ring cap N 11, PE, bluecenter	100 pz
LLG04008261	LLG-Snap ring cap N 11, PEtransparent,c	100 pz
LLG06073555	PE Snap Ring Seals ND11,ready assembled	100 pz
LLG06242212	LLG-Snap ring cap N 11blue, hole, s/100	100 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5ª Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

