



Septums LLG para tapones roscados ND13

LLG Labware

Septums LLG para tapones roscados ND13



Code	Description	Packaging
LLG04008263	LLG-Septa N 12, Silicone white/PTFE/100	100 pz
LLG07058143	LLG-Septa N 12, PTFE virginal, whit/100	100 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5ª Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

