



Viales Shell LLG con tapón de PE (cuello pulido a fuego) ND8, ND12, ND15

LLG Labware

Viales Shell LLG con tapón de PE (cuello pulido a fuego) ND8, ND12, ND15



Code	Description	Packaging
LLG04008205	LLG-Shell Vial N 8, 1ml O.D.: 8,2mm,/100	100 pz
LLG04008248	LLG-Shell Vial N 122ml with stopper/100	100 pz
LLG06226501	2 ml flat bottom glass, brown 12 mm/ 1000	1000 pz.
LLG06280950	LLG-Flat bottom glasses 4 ml/100	100 pz.
LLG07300174	LLG-sample vials N 8-1.2 W, colour/100	100 pz
LLG07612017	LLG Vials tipo Shell N 8 + N 12 Capacit	1000 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

