



Viales Shell LLG con tapón de PE (cuello pulido a fuego) ND8, ND12, ND15

LLG Labware

Viales Shell LLG con tapón de PE (cuello pulido a fuego) ND8, ND12, ND15



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG04008205 | LLG-Shell Vial N 8, 1ml O.D.: 8,2mm,/100 | 100 pz |
| LLG04008248 | LLG-Shell Vial N 122ml with stopper/100 | 100 pz |
| LLG06226501 | 2 ml flat bottom glass, brown 12 mm/ 1000 | 1000 pz. |
| LLG06280950 | LLG-Flat bottom glasses 4 ml/100 | 100 pz. |
| LLG07300174 | LLG-sample vials N 8-1.2 W, colour/100 | 100 pz |
| LLG07612017 | LLG Vials tipo Shell N 8 + N 12 Capacit | 1000 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

