



Cierres para encapsulado ND20 LLG, aluminio, ensamblados

LLG Labware

Cierres para encapsulado ND20 LLG, aluminio, ensamblados



Code	Description	Packaging
LLG04001549	LLG-Aluminium crimp cap N 20, silve/100	100 pz
LLG04008275	LLG-Aluminium crimp cap N 20, gold, cen	100 pz
LLG07060477	LLG-Aluminium crimp cap N 20, silve/100	100 pz
LLG09003430	LLG-Aluminium crimp cap N 20, silver, c	100 pz
LLG09003434	LLG-Alu crimp cap N 20, silverhole,/100	100 pz
LLG09003454	LLG-Aluminium crimp cap N 20, silve/100	100 pz
LLG09003460	LLG-Aluminium crimp cap N 20,silver/100	100 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

