



**LLG-Crimp Seals ND20, Aluminium, ready assembled**

LLG-Crimp Seals ND20, Aluminium, ready assembled



Code	Description	Packaging
LLG04001549	Crimp Seals ND20, Aluminium, ready assembled, Caps silver, center hole, Septa Butyl grey / PTFE grey, Hardness 50° shore A, Thickness 3 mm	100 pz
LLG04008275	Crimp Seals ND20, Aluminium, ready assembled, Caps gold, center hole, Septa Pharma-Fix- Septa Butyl/PTFE, Hardness 50° shore A, Thickness 3 mm	100 pz
LLG09003430	Crimp Seals ND20, Aluminium, ready assembled, Caps silver, center hole, Septa Pharma-Fix- Septa Butyl/PTFE, Hardness 50° shore A, Thickness 3 mm	100 pz
LLG09003434	Crimp Seals ND20, Aluminium, ready assembled, Caps silver, center hole, Septa Silicone blue / PTFE colourless, Hardness 45° shore A, Thickness 3 mm	100 pz
LLG09003454	Crimp Seals ND20, Aluminium, ready assembled, Caps silver, center hole, Septa Butyl red / PTFE grey, Hardness 50° shore A, Thickness 3 mm	100 pz
LLG09003460	Aluminium crimp cap N 20,silver/100	100 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

