



## Picnómetros

### Picnómetros

Conformes a ISO 2811, DIN 53 217, ASTM D1475.



Code	Description	Packaging
LLG06205559	Pyknometer PR 5 50ml made of stainles/1	1 pz.
LLG06223204	Pyknometer PR 10, 100 ml made of stai/1	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

