



Picnómetros

Picnómetros

Conformes a ISO 2811, DIN 53 217, ASTM D1475.



Code	Description	Packaging
LLG06205559	Pyknometer PR 5 50ml made of stainles/1	1 pz.
LLG06223204	Pyknometer PR 10, 100 ml made of stai/1	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

