



Tapones bi-metal para encapsulado ND20 LLG, ensamblados, magnéticos

LLG Labware

Tapones bi-metal para encapsulado ND20 LLG, ensamblados, magnéticos



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG06234541 | LLG-Bimetal crimp cap N 20, red/sil/100 | 100 pz |
| LLG09003457 | LLG-Bimetal crimp cap N 20, red/sil/100 | 100 pz |
| LLG09003458 | LLG-Bimetal crimp cap N 20, red/sil/100 | 100 pz |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

