



Tapones bi-metal para encapsulado ND20 LLG, ensamblados, magnéticos

LLG Labware

Tapones bi-metal para encapsulado ND20 LLG, ensamblados, magnéticos



Code	Description	Packaging
LLG06234541	LLG-Bimetal crimp cap N 20, red/sil/100	100 pz
LLG09003457	LLG-Bimetal crimp cap N 20, red/sil/100	100 pz
LLG09003458	LLG-Bimetal crimp cap N 20, red/sil/100	100 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5ª Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

