



Tapones de acero para encapsulado LLG ND20, ensamblados, magnéticos

LLG Labware

Tapones de acero para encapsulado LLG ND20, ensamblados, magnéticos



Code	Description	Packaging
LLG06229635	LLG-Magnetic crimp cap N 20, gold, /100	100 pz
LLG06902419	LLG-Steel crimp caps N 20 TB/oA-M, gold	100 pz
LLG07850010	LLG-Magnetic crimp cap N 20, gold, /100	100 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

