



Septums LLG para tapones para encapsulado ND20

LLG Labware

Septums LLG para tapones para encapsulado ND20



Code	Description	Packaging
LLG04001550	LLG-Septa N 20 Silicone white/Alumi/100	100 pz
LLG04008273	LLG-Septa N 20, Butyl light grey/PT/100	100 pz
LLG04008274	LLG-Septa N 20 Silicone blue/PTFE c/100	100 pz
LLG07050202	LLG-Septa N 20 Silicone beige/PTFE /100	100 pz
LLG07071063	LLG-Septa N 20 Butyl dark grey/cent/100	100 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

