



**LLG-Septa for Crimp Caps ND20**

**LLG Labware**

**LLG-Septa for Crimp Caps ND20**



Code	Description	Packaging
LLG04001550	Septa for Crimp Caps ND20, Septa Silicone white / Aluminium foil silver, Hardness 50° shore A, Thickness 3 mm	100 pz
LLG04008273	Septa for Crimp Caps ND20, Septa Butyl grey / PTFE grey, Hardness 50° shore A, Thickness 3 mm	100 pz
LLG04008274	Septa for Crimp Caps ND20, Septa Silicone blue transp./ PTFE transparent, Hardness 45° shore A, Thickness 3 mm	100 pz
LLG07050202	Septa for Crimp Caps ND20, Septa Silicone white / PTFE beige, Hardness 45° shore A, Thickness 3.2 mm	100 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

