



Tapones de PE transparentes ND20 LLG y septums adecuados

LLG Labware

Tapones de PE transparentes ND20 LLG y septums adecuados



Code	Description	Packaging
LLG06227768	LLG-Cap N 20, PE, for crimp neck vi/100	100 pz
LLG09003543	LLG-Closure 22 mm, PE, height 8,4m/ 1000	1000 pz
LLG09003544	LLG-Closure 22 mm, PE, height 8.4m/ 1000	1000 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

