



Cierres de rosca ND24, (Cierres Ultrabond y cierres ensamblados ND24), tapones de rosca de PP ND24 (vacíos) y septums ND22 LLG

LLG Labware

Cierres de rosca ND24, (Cierres Ultrabond y cierres ensamblados ND24), tapones de rosca de PP ND24 (vacíos) y septums ND22 LLG



Code	Description	Packaging
LLG04008295	LLG-Screw cap N 24, PP,white, cente/100	100 pz
LLG04008296	LLG-Screw-cap N 24, PP,white, close/100	100 pz
LLG06233339	Screw-cap PP, white with holesilic/ 1000	1000 pz.
LLG09003541	LLG-Screw cap PP, white mounted cl/ 1000	1000 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

