



Cierres de rosca ND24, (Cierres Ultrabond y cierres ensamblados ND24), tapones de rosca de PP ND24 (vacíos) y septums ND22 LLG

LLG Labware

Cierres de rosca ND24, (Cierres Ultrabond y cierres ensamblados ND24), tapones de rosca de PP ND24 (vacíos) y septums ND22 LLG



Code	Description	Packaging
LLG04008295	LLG-Screw cap N 24, PP,white, cente/100	100 pz
LLG04008296	LLG-Screw-cap N 24, PP,white, close/100	100 pz
LLG06233339	Screw-cap PP, white with holesilic/ 1000	1000 pz.
LLG09003541	LLG-Screw cap PP, white mounted cl/ 1000	1000 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

