



Tapones a presión LLG ND18 y ND22, LDPE

LLG Labware

Tapones a presión LLG ND18 y ND22, LDPE

Para el sellado de viales con aro de presión.



Code	Description	Packaging
LLG06051403	Snap Caps, PE, ND18 and ND22, For ND18	100 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

