



Tapones de rosca LLG para viales con rosca ND15, ND18

LLG Labware

Tapones de rosca LLG para viales con rosca ND15, ND18



Code	Description	Packaging
LLG04678391	LLG-Screw caps 15 mm, black PP, closed,	100 pz



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

