



Tapones de rosca LLG para viales con rosca ND15, ND18

LLG Labware

Tapones de rosca LLG para viales con rosca ND15, ND18



Code	Description	Packaging
LLG04678391	LLG-Screw caps 15 mm, black PP, closed,	100 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

