



Septos de sustitución LLG para tapones de rosca ND15, ND18

LLG Labware

Septos de sustitución LLG para tapones de rosca ND15, ND18



Code	Description	Packaging
LLG09003536	Septa 12mm, Butyl red/PTFEgrey, 55/ 1000	1000 pz.
LLG09003538	Septa 16mm, Butyl red/PTFEgrey, 55/ 1000	1000 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

