

**Microtubos Safe-Lock, PP****Eppendorf SE****Microtubos Safe-Lock, PP**

- Un pequeño gancho del cierre del tapón se engancha en el reborde del microtubo
- Autoclavables estando abiertos (121 °C, 20 min.)
- Se pueden abrir y cerrar fácilmente con una mano
- Con superficie esmerilada de rotulación
- El gancho del tapón impide que salte por accidente (por ej. durante la termostatización)
- Máxima estabilidad mecánica para la centrifugación
- Graduados
- Los recipientes de reacción Safe-lock también están disponibles en calidad Eppendorf Biopur®
- Los recipientes de reacción con cierre de seguridad de 0,5 ml y 1,5 ml de capacidad pueden ser centrifugados hasta 30 000 x g



Code	Description	Packaging
LLG09409320	SAFE-LOCK reaction vessels0.5 ml,co/500	500 pz.
LLG09409331	SAFE-LOCK reaction vessels1.5 ml,c/ 1000	1000 pz.
LLG09409332	SAFE-LOCK reaction vessels,cap. 1./ 1000	1000 pz.
LLG09409333	SAFE-LOCK reaction vessels,cap. 1./ 1000	1000 pz.
LLG09409334	SAFE-LOCK reaction vessels,cap. 1./ 1000	1000 pz.
LLG09409335	SAFE-LOCK reaction vessels,cap. 1./ 1000	1000 pz.
LLG09409336	SAFE-LOCK reaction vessels,cap. 1./ 1000	1000 pz.
LLG09409341	SAFE-LOCK reaction vessels2 ml,col/ 1000	1000 pz.
LLG09409342	SAFE-LOCK reaction vessels,cap. 2 / 1000	1000 pz.
LLG09409343	SAFE-LOCK reaction vessels,cap. 2 / 1000	1000 pz.
LLG09409344	SAFE-LOCK reaction vessels,cap. 2 / 1000	1000 pz.
LLG09409345	SAFE-LOCK reaction vessels,cap. 2 / 1000	1000 pz.
LLG09409346	SAFE-LOCK reaction vessels,cap. 2 / 1000	1000 pz.
LLG09409352	SAFE-LOCK reaction vessels,cap. 0.5/500	500 pz.
LLG09409353	SAFE-LOCK reaction vessels,cap. 0.5/500	500 pz.
LLG09409354	SAFE-LOCK reaction vessels,cap. 0.5/500	500 pz.
LLG09409355	SAFE-LOCK reaction vessels,cap. 0.5/500	500 pz.
LLG09409356	SAFE-LOCK reaction vessels,cap. 0.5/500	500 pz.
LLG09409357	Safe-Lock reaction vessel 0,5 mlass/500	500 pz.
LLG09409358	Safe-Lock reaction vessels 1,5 mla/ 1000	1000 pz.



Code	Description	Packaging
LLG09409359	Safe-Lock reaction vessels 2,0 mla/ 1000	1000 pz.



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

