

**Butirómetro para leche Original Gerber****Gerber****Butirómetro para leche Original Gerber**

Etiqueta esmerilada sobre el bulbo, cuello liso, sin tapón, para cantidades de leche de 10,73-11 ml.



Code	Description	Packaging
LLG06224904	Milchbutyrometer O-8% F Glatthals, 510m	1 pz.
LLG06225776	Milk butyrometer 0-10%:0.1	1 pz.
LLG06226307	Milk butyrometer 0-9% F	1 pz.
LLG06236033	Milk butyrometer 0-12%:0,1	1 pz.
LLG06236845	Milk butyrometer 0-14%:0,2	1 pz.
LLG06266488	Milk butyrometer 0-16%:0,2	1 pz.
LLG06803740	Milk butyrometer 0-7% F not calibrated,	1 pz.
LLG09112565	Milk/skimmed milk butyrometer0-4%:0,05	1 pz.
LLG09112588	Milk butyrometer 0-6%:0,1without stoppe	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

