



Skim Milk Butyrometer

Gerber

Skim Milk Butyrometer

With frosted label on bulb, without stopper.



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG09112564 | Skimmed milk butyrometer 0-1%:0.02 F, e | 1 pz |
| LLG09112566 | Skim Milk Butyrometer, Milk: 43.5ml \pm 0.5ml, closed bulb, plain neck, Fat 0 - 0.5 %, Graduation 0,02 | 1 pz |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

