



Butirometros para leche descremada

Gerber

Butirometros para leche descremada

Con etiqueta esmerilada sobre el bulbo, sin tapón.



Code	Description	Packaging
LLG09112564	Skimmed milk butyrometer 0-1%:0.02 F, e	1 pz.
LLG09112566	Skimmed milk butyrometer Siegfeld	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

