



Accesorios para Butirómetros



Code	Description	Packaging
LLG06202638	***no longer available*** ***replacemen	1 pz
LLG06224906	Rubber stopper double conical Ø 12/14x3	1 pz.
LLG06234037	Glass pin for milk butyrometer	1 pz.
LLG06262314	Glass beaker for butyrometers "butter'	1 pz.
LLG06264750	Weighing funnel, glass	1 pz.
LLG09112573	Rubber stopper GERBAL-Mfor volumetric m	1 pz.
LLG09112574	Key for GERBAL-M	1 pz.
LLG09112575	Rubber stoppers, large30/17x22mm with h	1 pz.
LLG09112576	Rubber stoppers, small15/10x13mm	1 pz.
LLG09112577	"Glass beaker for butyrometers""cream""	1 pz.
LLG09112578	"Glass beaker for butyrometers""cheese""	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*